
IntegraPost®

Ultra-passive Integrating Core Retainer

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Intended Use:
IntegraPost is intended to provide retention for the overlying core build-up and final restoration when insufficient coronal tooth structure remains.

US / UK
INSTRUCTIONS FOR USE

Indications for Use:
The Integra Post is indicated where less than half of the coronal tooth structure is remaining subsequent to the completion of endodontic therapy. The purpose of the post is to retain the core build-up material and reinforce the remaining coronal tooth structure.

Product Description:
Premier® IntegraPost® is an ultra-passive, ultra-secure integrating endodontic post made of titanium alloy. IntegraPost features a round Flow-Thru-Head to optimize composite core materials and a parallel passive design that creates an adhesive buffer zone. The IntegraPost system conserves tooth structure and minimizes stress.

SKU	Description
3001105	Size 3
3001110	Size 4
3001115	Size 4.5
3001125	Size 5
3001130	Size 6
3001301	Integra Instruments – IP Carrier
3001305	Integra Instruments – IP Probe

Integra Post, Refill (Package of 10)

that the instruments are sufficiently immersed Note: review special instructions for cleaning solutions to be used.

- An Ultrasonic Tank can be used. Use the processing time recommended by the manufacturer of the detergent.
- Remove posts from the cleaning solution and post rinse them intensively with low contaminated and deionized water.
- Thoroughly dry the instruments after cleaning.

Automated Cleaning

- Place the posts in a cassette or any other tray system suitable to ensure that there is no contact between the items and place it in the automated washer disinfecter.
- Start the program.
- Remove the posts from the automated washer disinfecter after end of the program.
- Allow post drying step in a clean place.

Sterilization:

- Inspect for corrosion, damage, defects, residue or soil. Return the device for re-cleaning if any soil or residue is still visible.
- Dispose of device where there are signs of corrosion, defects, or damages. Decommission and dispose of the device accordingly.
- Place cleaned and dried devices into FDA approved sterilization pouches, wraps, or suitable sterilization containers. Do not place dissimilar metals within the same sterilization wrap or pouch.
- Sterilize per instructions below

Minimum cycle times for dynamic-air-removal steam sterilization cycles (Pre-Vacuum)

EXPOSURE TEMPERATURE	EXPOSURE TIME	DRYING TIMES
270degF (132degC)	4 minutes	Minimum 30 minutes

Minimum cycle times for Gravity Steam Sterilization

EXPOSURE TEMPERATURE	EXPOSURE TIME	DRYING TIMES
250degF (121degC)	30 minutes	Minimum 30 minutes

Storage:

- Store in dry area.
- Do NOT use if product is damaged.

For Further Information:
To obtain IFU, visit www.premierdentalco.com or call Premier at 610-239-6000.

Symbols:

A full explanation of symbols used on Premier Dental packaging and literature is located at premierdentalco.com/symbol-glossary.

Symbols	Symbol Meaning
LOT	Batch number
REF	Catalog Number
US	Date and Country of Manufacture
Manufacturer	Manufacturer
i	See IFU on premierdentalco.com
X	Single Use, Do Not Re-Use
Rx ONLY	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician
NON-Sterile	Non-Sterile
DO NOT USE	Do not use if package is damaged

Ordering Information:

CompCore™ AF Composite Core Material	
SyringeMix Stack: Dual-Cure SyringeMix contains two 9g auto-mix syringes, 15 mix-tips, 15 standard & 10 extra-fine intraoral tips	Shade White 3001512
	Shade A3 3001513
	Shade Blue 3001514
Cartridge Stack: Dual-Cure Intro Kit contains one 50g cartridge, type-25 gun, 4ml each IntegraBond & Auto-cure Activator, 25 mix-tips, 25 intraoral tip	Shade White 3001430
	Shade A3 3001431
Cartridge Refills Refill contains one 50g cartridge, 10 ACCOR® matrices, 25 mix-tips, 25 intraoral tips	Shade White 3001432
	Shade A3 30014313
CompCore Auto-Cure/Tub Hand-Mix Tub contains catalyst & base, mixing pads, spatulas	Shade Natural (28g) 3001405
	Shade Natural (50g) 3001404

Accessories

SyringeMix Tips & Intraoral Tips (35)	3001418
SyringeMix XXF Extra-fine Intraoral Tips (50)	3001468
Applicator Gun, Type-25	3001416
Cartridge Mix-Tips & Intra-oral Tips (35)	3001417

Integra Adhesives

IntegraBond™ Adhesive Kit	3001421
Premier® Etch 37%	3001422
IntegraBond 7ml	3001411
IntegraBond 4ml	3001413
Auto-Cure Activator 7ml	3001412

IntegraPost

IntegraPost Refill (10)	Size 3 3001105
	Size 4 3001110
	Size 4.5 3001115
	Size 5 3001125
	Size 6 3001130

Cure-Thru IntegraPost Refill (10)	Size 4 3001610
	Size 4.5 3001611
	Size 5 3001612
	Size 6 3001613

CLEANING & STERILIZATION

Clean and sterilize, prior to use. Do NOT re-process.

Cleaning:

Manual Cleaning

- Submerge the post in a cleaning solution. Soak for the recommended soaking time in the cleaning solution and make sure

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